



### 1100 Fact!

The 1100 can be installed quickly – our dealers are trained to be efficient; stairlift installation should be quick and pain-free and shouldn't interfere with your day.



Vanilla Beige

## A Stairlift Suited for Any Interior

We designed the 1100 to be a smart, stylish complement to everything else in your home. The contemporary design is welcome in any home décor and is available in a warm vanilla beige upholstery.

Featuring our comfortable Style Seat in a leatherette finish, the seat surface is easily wipeable and cleans effortlessly.

The ultra slim rail beautifies the 1100, simply there is no other rail on the market that competes with its overall look and sleek profile, a non-obtrusive design that quietly blends in with your staircase.

# Frequently Asked Questions

We strive for 100% customer satisfaction. That is why we go the extra mile to offer you the highest standards of service and answer your questions in the best way possible. We partner with the industry's leading dealers that have a track record of exceptional local support.

We have listed some of our customers' most frequently asked questions. Should you have further questions, please contact your local dealer for a consultation.



## Can a stairlift be installed on my staircase?

Handicare offers a wide selection of stairlifts for straight and curved stairs. With a variety of lift models we are able to provide a suitable solution for each setting regardless of whether the stairs are narrow, steep, curved or long, indoor or outdoor. Handicare always has a solution for your staircase.

## Does it take major renovation to install a stairlift in my home?

No! The 1100 Straight Stairlift is simply installed onto the stair treads. Because it is fixed to the staircase, you don't have to worry about any further modification to your home.

## What happens to my stairlift during a power outage?

A power failure is no reason to worry. Your Handicare stairlift will use its batteries to remain operational and can still be used for at least 10 rides up and down the stairs. The 1100 features Constant Charge, providing charging anywhere on the track so it is ready for use 100% of the time.

## How do I know when batteries need to be replaced?

The 1100 is equipped with a fault display behind the power pack, simply review the code and call your local dealer who can easily look up the code and inform you if the battery needs to be replaced.

Generally the batteries will last for three years depending on usage.

## How do I enter and exit from a stairlift?





Entering and exiting the 1100 is made easy, simply unfold the seat and footrest, sit down while securing the seat belt. When it's time to exit your 1100 you can use easily swivel the seat to face your landing for a controlled departure.

## Can others still use my stairs?

Compact and ultra slim, the 1100 can easily be folded up when not in use to allow for easy pass-through for family members and friends.

# 1100 Straight Stairlift

## Specifications



### 1100 Style Seat Dimensions

#### Safe Working Load (SWL)

30–45°: 309 lb (140 kg)

46–50°: 287 lb (130 kg)

#### Staircase Angle

30–50°

#### Speed

0.5 ft / sec (0.15 m / sec)

Some of the dimensions depend on the characteristics of the staircase in the home. Handicare reserves the right to alter designs and specifications without prior notice.

<b>A</b> Swivel radius, minimum	25.4 in (645 mm)
<b>B</b> Wall to inner seat back	2.4 in (60 mm)
<b>C</b> Seat depth	15.4 in (390 mm)
<b>D</b> Folded width to footplate, minimum	12.6 in (320 mm)
<b>E</b> Open width to edge footplate, minimum	20.9 in (530 mm)
<b>F</b> Armrest width, external	22.6 in (575 mm)
<b>G</b> Armrest width, internal	17.7 in (450 mm)
<b>H</b> Seat back height	15.7 in (400 mm)
<b>I</b> Footplate to seat height	16.9–20.5 in (430–520 mm)
<b>J</b> Footplate height, minimum	2.75 in (70 mm)
<b>K</b> Track intrusion into staircase, minimum	6.3 in (160 mm)

Are you ready to take back your independence with the 1100?


**STAIRLIFTS**  
**MADE SIMPLE**  
**951-526-2073**